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# OPERATING INSTRUCTIONS

## I. APPLICATION

The superhigh-pressure quartz spherical mercury lamps, type ДРШ, are powerful concentrated sources of radiation in the visible and ultra-violet spectrum range.

The lamps are intended for use in optical equipments and for laboratory research.

## II. SPECIFICATIONS

The lamps, type ДРШ, are intended for operation on 127 V A.C. (lamps ДРШ-250-3) and on 220 V A.C. (lamps ДРШ-500М and ДРШ-1000) in series with chokes. During the first 10 to 15 minutes after switching-on electrical parameters of the lamps vary (unstable period) and then become constant at supply voltage.

Electrical parameters of the lamps operating in an A.C. circuit and their dimensions are given in Table 1 and Fig. 1.

Table 1

Type of lamp	Wattage, W	Lamp voltage, V			Luminous flux, lm		Brightness in beam centre, M-nit	Dimensions, mm	
		rated	max.	min.	rated	min.		L <sub>max</sub>	b <sub>max</sub>
ДРШ-250-3	250	72	87	57	12 500	10 000	100	145	40
ДРШ-500М	500	75	92	60	22 100	18 000	130	190	45
ДРШ-1000	1000	90	105	75	53 000	42 500	120	232	58

### III. CONSTRUCTION AND PRINCIPLE OF OPERATION

The lamp, type ДРШ, is a spherical bulb made of quartz glass, with two main electrodes sealed into the bulb. The electrodes are arranged diametrically opposite to each other at a distance of 3.9 to 4.5 mm. An auxiliary electrode sealed into the bulb side facilitates lamp ignition. The bulb is filled with a strictly accurate amount of mercury and argon. An arc discharge occurs in mercury vapour, producing a high degree of brightness. A schematic diagram for connecting the lamp to an A.C. circuit is shown in Fig. 2.

The lamps are started by a magneto (with a spark length of 15 to 20 mm) which is engaged by push-button switch  $K$  for a period of not over 1 minute. After an arc discharge takes place between the current-carrying electrodes, the magneto is disengaged.

A lead from the magneto should be connected to the auxiliary electrode base.

#### Choke Characteristics

Table 2

Type of lamp	Operating conditions		Starting conditions	
	choke voltage, $V$	current, A	choke voltage, $V$	current, A
ДРШ-250-3	93	$40 \pm 0.2$	1.7	$6.5 \div 8.0$
ДРШ-500М	190	$7.5 \pm 0.3$	220	$9.0 \div 11.0$
ДРШ-1000	190	$12 \pm 0.4$	220	$15 \div 17$

### IV. SAFETY RULES

The fitting body should be grounded, and electrical connections should ensure reliable contact.

The bulb is attached by lugs or bases of current-carrying electrodes.

The bulb can operate in a restricted space (enclosure) provided the size of the enclosure and its ventilation are that the air temperature measured at a distance of 6 cm from the sides, does not exceed +250°C (steady conditions). An interval required in-between repeated switchings does not exceed 6 minutes.

The quartz glass of the bulb readily loses its transparency if dirty. Therefore, it is necessary to wipe the bulb with a piece of cotton-wool soaked in alcohol each time it was touched by hands or soiled with dust.

Personnel, equipment and inflammable materials in the vicinity of the lamp should be protected against hot bulb fragments in case of a burst of the lamp.

Measures should be taken to protect attending personnel against powerful ultra-violet radiation of the lamp.

The normal operating position of the lamp is vertical with the auxiliary electrode facing upward (the maximum permissible deviation from the vertical being  $10^\circ$ ). The lamp starts and operates reliably at ambient air temperatures not below  $+15^\circ\text{C}$  and supply voltage not below 115 V for the lamps of type ДРШ-250-3, 200 V for the lamps ДРШ-500M, and 205 V for the lamps ДРШ-1000.