


MICROSCOPES



MASHPRIBORINTORG
SSSR · MOSKVA





The Soviet optico-mechanical industry turns out a wide range of different microscope types which are used with great success at factories, educational establishments, and scientific centres both in the Soviet Union and abroad.

Using the latest advances in optical instrument-making, a large contingent of research workers, engineers, and technicians are continually perfecting existing microscopes, improving their optical systems and components, and developing fundamentally new designs and types of microscopes for a great variety of applications.

Recently a number of microscopes for spectral and other microscopic investigations have been developed (МСЛ-2, МЛК-1, МСК-1, МИК-3, etc.). At the same time, some microscopes have been brought up to date (among them the "Biolum" series of microscopes, ММУ-3, МИМ-9, МИК-4, ММР-2Р, etc.), which has made it possible not only to improve quality and increase productivity, but also to extend the service life of the instruments.

A wide choice of all kinds of objectives and eyepieces along with the high quality of Soviet-made optics known throughout the world, high technical characteristics, convenience and reliability in operation, and functional styling are all features of the Soviet-made microscopes designed to handle various problems in science and technology.

Please send all inquiries for additional information about Soviet-made instruments and orders for them to:

121200 Moscow, Smolenskaya-Sennaya, 32-34,
V/O MASHPRIBORINTORG

Telegrams: Moscow Mashpribor

Telephone: 244-27-75

Telex: 7235, 7236

CONTENTS

BIOLOGICAL MICROSCOPES

	Page
Type ШМ-1 classroom microscope	4
Type МБУ-4 simplified biological microscope	5
Type МБУ-6 trichina microscope	6
"BIOLAM P-1" practical biological microscope	7
"BIOLAM P-2" practical biological microscope	8
"BIOLAM P-3" practical biological microscope	9
"BIOLAM P-4" practical biological microscope	10
"BIOLAM P-5" practical biological microscope	11
"BIOLAM P-6" practical biological microscope	12
"BIOLAM C-1" students' biological microscope	13
"BIOLAM C-2" students' biological microscope	14
"BIOLAM C-3" students' biological microscope	15
"BIOLAM Д-1" travelling biological microscope	16
"BIOLAM Д-2" travelling biological microscope	17
Type МБИ-3 biological research microscope	18
Type МБИ-6 universal biological microscope	20
Type МББ-1А large biological microscope	22
Type МБИ-11 universal research microscope	24
Type МБИ-13 inverted biological microscope	26
Type МБИ-15 universal research microscope	28
"BIOLAM Д-3" travelling biological microscope	30

STEREOSCOPIC MICROSCOPES

Type МБС-2 stereoscopic microscope	32
Type МПС-2 stereoscopic polarizing microscope	34
Type БМ-51-2 binocular microscope	36
Type МБС-200 stereoscopic microscope	37
Type МССО stereoscopic microscope	38

METALLOGRAPHIC MICROSCOPES

Type ММУ-3 simplified metallographic microscope	40
Type ММР-2 metallographic shop microscope	41
Type ММР-2Р metallographic microscope	42
Type МИМ-9 universal metallographic microscope	44

POLARIZING MICROSCOPES

Type МИН-8 research polarizing microscope	48
Type МПД-1 travelling polarizing microscope	50
Type МИН-9 ore microscope	52

FLUORESCENCE MICROSCOPES

Type МЛ-2А fluorescence microscope	54
Type МЛ-2Б fluorescence microscope	56
Type МЛ-2В fluorescence microscope	57
Type МЛ-3 fluorescence microscope	58

Type МЛД-1 travelling fluorescence microscope	60
Type МЛК-1 contact fluorescence microscope	62
Type МЛ-4 fluorescence microscope	64

INFRA-RED AND ULTRA-VIOLET MICROSCOPES

Type МИК-4 infra-red microscope	66
Type МУФ-6 ultra-violet microscope	68
Type МИК-3 infra-red stereoscopic microscope	70

MICROPROJECTION, PHOTOMICROGRAPHY AND CINE-MICROGRAPHY OUTFITS

Type ТМП trichina microprojector	72
Type ФМН-2 photo macro- and micrography outfit	74
Type ФМН-3 photomicrographic outfit	76
Type МКУ-2 cinemicrography outfit	77

MICROSCOPES FOR ADVANCED MICROSCOPIC INVESTIGATIONS

Type МСЛ-2 microspectral analysis outfit	80
Type МСК-1 criminal investigation comparison microscope	82
Type МБИН-4 interference microscope	84

MEASURING MICROSCOPES

Type ММИ-2 small tool-room microscope	86
Type БМИ-1 large tool-room microscope	88
Type УИМ-23 universal measuring microscope	90
Type УИМ-24 universal measuring microscope	91
Type МИР-2 reading microscope	92
Type МПБ-2 reading microscope	93
Type МПВ-1 hardness-test microscope	94
Type МИР-12 microscope for spectrogram measurements	95

SURFACE INSPECTION INSTRUMENTS

Type МИИ-10 microinterferometer	98
Type МИИ-11 multi-beam microinterferometer	100
Type МИИ-4 interference microscope	102
Type МИИ-12 micromprofilometer	103

Biological Microscopes



CLASSROOM MICROSCOPE WM-1

The WM-1 class-room microscope is the simplest biological microscope designed for bright-field, transmitted-light work at low magnifications.


The WM-1 class-room microscope is an aid in the study of natural sciences at secondary and technical schools and other educational establishments and is intended to show a variety of preparations to students and to examine minute details of objects invisible to the unaided eye.

In addition to direct observation in the microscope, it may be supplemented with a variety of attachments or simple fixtures to project microscopic objects onto a viewing screen, to take their photographs on photographic plates (photomicrography), to draw sketches, or to measure objects, etc.

The microscope comes complete with two achromatic objectives (8×0.20 and 20×0.40) and three Huygens eyepieces (7×, 10× and 15×).

BASIC CHARACTERISTICS

Total magnification of microscope	56 to 300×
Angle of tilt of microscope	≤45°
Overall dimensions, mm	130×180×270
Mass, kg	2.4



4

SIMPLIFIED BIOLOGICAL MICROSCOPE MBY-4

The MBY-4 is a simplified biological microscope mainly intended for transmitted light, bright-field work and also for the examination of opaque objects at low magnifications. The MBY-4 microscope is designed for use at general and technical secondary schools and other educational establishments as an aid in the study of natural sciences, at botanical, agricultural, zoological and medical laboratory, and also at factory laboratories.

The microscope comes complete with two achromatic objectives (8×0.20 and 20×0.40) and three Huygens eyepieces (7×, 10× and 15×).

BASIC CHARACTERISTICS

Magnification of microscope	56 to 300×
Drawtube travel per revolution of micrometer drum, mm . . .	0.1
Drawtube travel per division of micrometer drum rotation, mm	0.002
Overall dimensions, mm	130×205×305
Mass, kg	2.7



5

TRICHINA MICROSCOPE МБУ-6

The МБУ-6 trichina microscope is a simplified biological microscope intended for the detection and examination of trichina and other parasites in meat and meat products.


The microscope is designed for transmitted-light, dark-field work.

The microscope is widely used by meat-packing laboratories, livestock breeders, etc.

The microscope comes complete with two achromatic objectives (3.7×0.11 and 8×0.20) and one Huygens eyepiece (7×).

BASIC CHARACTERISTICS

Magnification of microscope	26× and 56×
Barrel travel per control knob revolution, mm	20
Stage dimensions, mm	100×170
Overall dimensions, mm	190×220×325
Mass, kg	2.86



6

PRACTICAL BIOLOGICAL MICROSCOPE "BIOLAM P-1"

The "Biolam P-1" practical microscope is intended for the transmitted-light, bright-field examination of transparent objects for educational and laboratory purposes in medicine, biology, zoology, and allied sciences.

With the "Biolam P-1" microscope, objects are examined through a monocular attachment.

The light system is a demountable mirror which is used for work with both natural (daylight) illumination and with a type ОН-19 illuminator supplied on separate order.

The microscope comes complete with three achromatic objectives (8×0.20, 40×0.65, and 90×1.25) and two compensating eyepieces (7× and 15×).

The "Biolam P-1" microscope is widely used by secondary and higher educational establishments, research institutions, and medical laboratories.

BASIC CHARACTERISTICS

Magnification of microscope	56 to 1350×
Focusing range, mm:	
coarse focusing	40
fine focusing	2
Overall dimensions, mm	230×140×350
Mass, kg	3.1



7

PRACTICAL BIOLOGICAL MICROSCOPE "BIOLAM P-2"

The "Biolam P-2" practical biological microscope is intended for the transmitted-light, bright-field examination of transparent objects and is widely used by educational and research establishments in studies in the fields of medicine, biology, zoology, and allied sciences.


With the "Biolam P-2" microscope, objects are examined through a monocular attachment.

The illuminating system of the microscope incorporates a type ON-32 simplified illuminator fitted into a locating socket in the base of the microscope, and a demountable mirror for work with natural (daylight) illumination or with a type ON-19 illuminator supplied on separate order.

The microscope comes complete with four achromatic objectives (8×0.20 , 20×0.40 , 40×0.75 , and 90×1.25) and two compensating eyepieces ($7\times$ and $15\times$).

BASIC CHARACTERISTICS

Magnification of microscope	56 to 1350 \times
Focusing range, mm:	
coarse focusing	40
fine focusing	2
Overall dimensions, mm	230 \times 140 \times 350
Mass, kg	3.1



8

PRACTICAL BIOLOGICAL MICROSCOPE "BIOLAM P-3"

The "Biolam P-3" practical biological microscope is intended for the transmitted-light, bright-field examination of transparent objects in investigations in the fields of medicine, biology, zoology, and allied sciences at educational and research establishments, in medical laboratories, etc.

The optical system of the microscope has provisions for observation of objects through a type AY-12 binocular attachment.

The illuminating system is a demountable mirror which can be used with both natural (daylight) illumination and with a type ON-19 illuminator supplied on separate order.

The microscope comes complete with three achromatic objectives (8×0.20 , 40×0.65 , and 90×1.25) and two paired compensating eyepieces ($7\times$ and $10\times$).

BASIC CHARACTERISTICS

Magnification of microscope	84 to 1350 \times
Focusing range, mm:	
coarse focusing	40
fine focusing	2
Overall dimensions, mm	230 \times 140 \times 350
Mass, kg	3.7



9

PRACTICAL BIOLOGICAL MICROSCOPE "BIOLAM P-4"

The "Biolam P-4" practical biological microscope is designed for the transmitted-light, bright-field examination of transparent objects.

The optical system of the microscope has provisions for observations of microscopic objects through a type AY-12 binocular attachment.


The illuminating system of the microscope incorporates a type ON-35 special-purpose demountable illuminator fitted into a locating socket in the microscope base, and a demountable mirror for work with both natural (daylight) illumination and with a type ON-19 illuminator supplied on separate order.

The microscope comes complete with three achromatic objectives (20×0.40 , 40×0.65 , and 90×1.25), one planachromatic objective (9×0.20), two paired compensating eyepieces ($7\times$ and $10\times$), and one measuring compensating eyepiece with a magnification of $7\times$.

The "Biolam P-4" microscope is used by educational and research establishments, medical laboratories in studies in the fields of biology, zoology and allied sciences.

BASIC CHARACTERISTICS

Magnification of microscope	94.5 to 1350 \times
Focusing range, mm:	
coarse focusing	40
fine focusing	2
Overall dimensions, mm	230 \times 350 \times 350
Mass, kg	4.5



PRACTICAL BIOLOGICAL MICROSCOPE "BIOLAM P-5"

The "Biolam P-5" practical biological microscope is designed for the transmitted-light, bright-field examination of transparent objects in various investigations in the fields of medicine, biology, zoology and allied sciences.


With the "Biolam P-5" microscope, objects are examined through a type AY-12 binocular attachment. The illuminating system is a demountable mirror which can be used with both natural (daylight) illumination and a type ON-19 illuminator supplied on separate order.

The microscope comes complete with three achromatic objectives (8×0.20 , 40×0.65 and 90×1.25), and two paired compensating eyepieces ($7\times$ and $10\times$).

The "Biolam P-5" microscope is widely used by educational and research establishments and medical laboratories.

BASIC CHARACTERISTICS

Magnification of microscope	84 to 1350 \times
Focusing range, mm:	
coarse focusing	40
fine focusing	2
Overall dimensions, mm	230 \times 140 \times 350
Mass, kg	4.5



PRACTICAL BIOLOGICAL MICROSCOPE

"BIOLAM P-6"

The "Biolam P-6" practical biological microscope is designed for the transmitted-light, bright-field examination of transparent objects and is widely used by educational and research establishments and medical laboratories in investigations in the fields of biology, zoology, medicine and allied sciences.


The optical system of the microscope has provisions for observation of microscopic objects through a type AY-12 binocular attachment.

The illuminating systems of the microscope incorporates a type ON-35 special-purpose demountable illuminator fitted into a locating socket in the microscope base. Besides, the standard equipment of the microscope includes a demountable mirror for work with natural (daylight) illumination or with a type ON-19 illuminator available on separate order.

The microscope comes complete with three achromatic objectives (20×0.40 , 40×0.65 , and 90×1.25), two paired compensating eyepieces ($7\times$ and $10\times$), and one measuring compensating eyepiece with a magnification of $7\times$.

BASIC CHARACTERISTICS

Magnification of microscope	84 to 1350 \times
Focusing range, mm:	
coarse focusing	40
fine focusing	2
Overall dimensions, mm	230 \times 140 \times 350
Mass, kg	4.5



12

STUDENTS' BIOLOGICAL MICROSCOPE

"BIOLAM C-1"

The "Biolam C-1" students' biological microscope is designed for the transmitted-light, bright-field examination of transparent objects in investigations in the fields of medicine, zoology, biology, and allied sciences.

With the "Biolam C-1" microscope, objects are examined through a monocular attachment.


The illuminating system is a demountable mirror which can be used with both natural (daylight) illumination and a type ON-19 illuminator supplied on separate order.

The microscope comes complete with three achromatic objectives (8×0.20 , 40×0.65 , and 90×1.25), and two compensating eyepieces ($7\times$ and $15\times$).

The "Biolam C-1" microscope is used by a secondary and higher educational establishments, research institutions, medical and other laboratories.

BASIC CHARACTERISTICS

Magnification of microscope	56 to 1350 \times
Focusing range, mm:	
coarse focusing	40
fine focusing	2
Overall dimensions, mm	230 \times 140 \times 350
Mass, kg	2.9



13

STUDENTS' BIOLOGICAL MICROSCOPE

"BIOLAM C-2"

The "Biolam C-2" students' biological microscope is designed for the transmitted-light, bright-field examination of transparent objects in investigations in the fields of medicine, zoology, biology and allied sciences.


With the "Biolam C-2" microscope, objects are viewed through a monocular attachment.

The illuminating system of the microscope incorporates a type OM-32 simplified demountable illuminator fitted into a locating socket in the microscope base. Besides, the standard equipment of the microscope includes a demountable mirror for work with both natural (daylight) illumination and a type OM-19 illuminator supplied on separate order.

The microscope comes complete with three achromatic objectives (8×0.20 , 20×0.40 , and 40×0.65), one planachromatic objective (3.5×0.10), two compensating eyepieces ($10\times$ and $15\times$), and one Huygens eyepiece with a magnification of $7\times$.

BASIC CHARACTERISTICS

Magnification of microscope	24.5 to 600×
Focusing range, mm:	
fine focusing	2
coarse focusing	40
Overall dimensions, mm	230×140×350
Mass, kg	2.9



14

STUDENTS' BIOLOGICAL MICROSCOPE

"BIOLAM C-3"

The "Biolam C-3" students' biological microscope is designed for the transmitted-light, bright-field examination of transparent objects in investigations in the fields of medicine, biology, zoology, and allied sciences.

The optical system of the microscope has provisions for observation of objects through a type AY-12 binocular attachment.

The illuminating system is a demountable mirror which can be used for work with both natural (daylight) illumination and a type OM-19 illuminator supplied on separate order.

The microscope comes complete with four achromatic objectives (8×0.20 , 40×0.65 , 85×1.0 and 90×1.25), and two paired compensating eyepieces ($7\times$ and $10\times$).

BASIC CHARACTERISTICS

Magnification of microscope	84 to 1350×
Focusing range, mm:	
coarse focusing	40
fine focusing	2
Overall dimensions, mm	230×140×350
Mass, kg	3.45



15

TRAVELLING BIOLOGICAL MICROSCOPE

"BIOLAM Д-1"

The "Biolam Д-1" travelling biological microscope is designed for the transmitted-light, bright-field examination of transparent objects in educational and research work in the fields of medicine, biology, zoology, and allied sciences under expeditionary conditions.

With the "Biolam Д-1" microscope, objects are viewed through monocular attachments. The standard equipment of the microscope includes a demountable mirror for work with both natural (daylight) illumination and a type ОИ-19 illuminators supplied on separate order.

The microscope comes complete with three achromatic objectives (40×0.75, 40×0.65, and 90×1.25), one planachromatic objective (9×0.20), two compensating eyepieces (10× and 15×), and one Huygens eyepiece with a magnification of 7×.

BASIC CHARACTERISTICS

Magnification of microscope	63 to 1350×
Focusing range, mm:	
coarse focusing	40
fine focusing	2
Overall dimensions, mm	230×140×350
Mass, kg	2.9



TRAVELLING BIOLOGICAL MICROSCOPE

"BIOLAM Д-2"

The "Biolam Д-2" travelling biological microscope is designed for the transmitted-light, bright-field examination of transparent objects in educational and research work in the fields of medicine, zoology, biology and allied sciences under expeditionary conditions.


With the "Biolam Д-2" microscope, objects are viewed through a monocular attachment.

The illuminating system of the microscope employs a type ОИ-32 simplified demountable illuminator fitted into a locating socket in the microscope base. Besides, the standard equipment of the microscope includes a demountable mirror for work with natural (daylight) illumination and type ОИ-19 illuminator supplied on separate order.

The microscope comes complete with three achromatic objectives (8×0.20, 40×0.65, and 90×1.25), two compensating eyepieces (10× and 15×), and one Huygens eyepiece with a magnification of 7×.

BASIC CHARACTERISTICS

Magnification of microscope	56 to 1350×
Focusing range, mm:	
coarse focusing	40
fine focusing	2
Overall dimensions, mm	230×140×350
Mass, kg	2.9



BIOLOGICAL RESEARCH MICROSCOPE

МБИ-3

The МБИ-3 biological research microscope is the simplest type of research microscope designed for the transmitted-light, bright- and dark-field examination of transparent objects with both direct and oblique illumination. The microscope is employed by medical and biological laboratories in the investigations of histological sections, microbiological preparations, to count formed blood elements, and also in other fields of science, technology and industry.


The microscope is convenient to use owing to an inclined binocular attachment which has a magnification of $1.5\times$ and a horizontally arranged stage. The binocular attachment may be replaced with a straight monocular drawtube to take photographs of the objects viewed.

All work on the microscope can be carried out with both natural (daylight) and artificial illumination. In the latter case, it will be good practice to employ types ОИ-19 and ОИ-31 illuminators which are not included in the standard equipment of the microscope, but may be had on separate order.

The microscope comes complete with four apochromatic objectives (10×0.30 , 20×0.65 , 60×0.7 to 1.0 , and 90×1.30), one achromatic objective (90×1.25), and six sizes of compensating eyepieces ($5\times$, $7\times$, $10\times$, $15\times$, $20\times$).

BASIC CHARACTERISTICS

Magnification of microscope	75 to $1350\times$
Least division of micrometer drum, mm	0.002
Overall dimensions, mm	$130\times 220\times 350$
Mass, kg	5.1



UNIVERSAL RESEARCH MICROSCOPE МБИ-6

The МБИ-6 microscope is a universal research microscope designed for work under laboratory conditions at medical, biological, physical, and chemical research institutions and may also be used in many fields of technology.

The МБИ-6 microscope can be used for examination and photography of objects by transmitted light with direct and oblique illumination; for dark-field, polarized-light, phase-contrast, and also for reflected-light bright- and dark-field work.

For photography work, the microscope can take attachable photographic cameras, namely a 35-mm film camera with 24×36 mm negatives, or a 9×12 cm plate camera.


The microscope comes complete with four planachromatic objectives for transmitted-light work (3.5×0.10, 9×0.20, 20×0.40, and 40×0.65), three apochromatic objectives (10×0.65, 60×0.70 to 1.0 and 90×1.30), five achromatic objectives for phase-contrast work (10×0.30, 20×0.40, 40×0.65, 90×1.25, and 40×0.75), one apochromatic objective for phase-contrast work (70×1.23), and four achromatic epi-objectives for reflected-light work (9×0.20, 20×0.40, 40×0.65, and 95×1.0).

The standard equipment of the microscope also includes five sizes of compensating eyepieces (5×, 7×, 10×, 15×, and 20×), two Huygens eyepieces (7× and 10×), and five photographic eyepieces (7×, 10×, 12.5×, 15×, and 20×).

BASIC CHARACTERISTICS

Magnification of microscope:

with transmitted light	17.5 to 2250×
with reflected light	45 to 2375×
Least increment in object displacement, mm	0.1
Maximum height of object in reflected-light work, mm	115
Overall dimensions, mm	1100×750×1440
Mass, kg	100



LARGE BIOLOGICAL MICROSCOPE MББ-1A

The MББ-1A large biological microscope is designed for the examination of transparent objects by transmitted light against bright and dark field with direct and oblique illumination.

The MББ-1A microscope is used by medical and biological research laboratories in the investigations of hystological sections, microbiological preparations and similar objects and is especially recommended for work which involves living low-contrast objects.

The contrast of unstained preparations may be improved and objects in vivo may be investigated by the phase contrast method. As an alternative, investigations may be carried out with polarized light.

For observations by reflected light against a dark and bright field, it is recommended to employ a type ОИ-21 illuminator which is not regularly supplied with the outfit but may be had on separate order.

The microscope comes complete with eight objectives as follows: one planachromatic objective (3.5×0.10), six apochromatic objectives (10×0.30 ; 20×0.65 ; 40×0.95 ; $60\times 0.7-1.0$; 70×1.23 ; and 90×1.30), one achromatic objective (40×0.75) and two sizes of compensating eyepieces.

The microscope is fitted with a type АУ-26 demountable binocular attachment with a triple revolving nosepiece which makes it possible to change the magnification of the attachment in steps ($1.1\times$, $1.6\times$, and $2.5\times$).

Photographs of microscopic objects can be taken with type МФН photomicrographic attachments supplied on separate order.

BASIC CHARACTERISTICS

Magnification of microscope	27 to 2250 \times
Stage:	
range of rotation	0-180°
range of vertical displacement, mm	0-50
range of transverse displacement, mm	0-40
forward and backward displacement	
range of mechanical stage, mm	0-60
Accuracy of stage and mechanical stage	
displacement readings, mm	0.1
Overall dimensions, mm	200 \times 450 \times 400
Mass, kg	6.5



UNIVERSAL RESEARCH MICROSCOPE МБИ-11

The МБИ-11 is a universal research microscope used by the laboratories of medical, biological, zoological and other research institutions. The microscope is adapted for the observation of objects by transmitted light against a bright field with direct and oblique illumination, against a dark field by the phase contrast method, by reflected light against a dark or bright field, and also with mixed illumination, that is, with the object simultaneously illuminated from below, through the condenser, and from above, through the objective.

The microscope is fitted with a type АУ-26 demountable binocular attachment whose magnification can be adjusted to 1.1 \times , 1.6 \times , and 5 \times . Instead of the binocular attachment, the microscope can be fitted with a type МФН photomicrographic attachment to take photographs of microscopic objects.

The microscope comes complete with six apochromatic objectives for transmitted-light work (10 \times 0.30; 20 \times 0.65; 40 \times 0.95; 60 \times 0.7-1.0; 90 \times 1.23; 90 \times 1.30), and four epi-objectives for reflected-light work (10 \times 0.20; 21 \times 0.40; 40 \times 0.65; and 95 \times 1.0). The microscope is also supplied with four sizes of compensating eyepieces (7 \times , and 10 \times), and a Huygens eyepiece with a magnification of 7 \times .

BASIC CHARACTERISTICS

Magnification of optical train of microscope:

by transmitted light (at maximum numerical aperture of 1.30)	92.4 to 2700 \times
by reflected light	69 to 2375 \times
stage displacement range in two mutually perpendicular directions, mm	55 \times 40
accuracy of stage displacement readings, mm	0.1
maximum working distance, mm	90
overall dimensions, mm	470 \times 200 \times 530
mass, kg	26.2



INVERTED BIOLOGICAL MICROSCOPE МБИ-13

The МБИ-13 inverted biological microscope is intended for visual observations, taking photographs and motionpicture films of living biological objects enclosed in a culture medium at a fixed temperature of $+37^{\circ}\text{C}$ or at any other predetermined temperature in the range from $+25$ to $+42^{\circ}\text{C}$.

The microscope is adapted for observations by transmitted light against a bright or dark field, with polarized light, and by the phase contrast method.

The instrument is a stationary outfit consisting of an inverted biological microscope proper, a constant temperature cabinet, and a picture-taking compartment containing single-shot and continuous-sequence mechanisms.

Photographs are taken with a ZORKY camera, and continuous sequences with type KCP-1M or 16CP cine cameras which can be actuated by both their own drives and the single-shot and continuous-sequence mechanisms.

The microscope can be used for investigations in oncology, virology, experimental histology, cytology, and other fields of medicine.

The microscope comes complete with two planachromatic objectives (3.5×0.10 and 9×0.20), three achromatic objectives (20×0.40 ; 40×0.65 ; and 60×0.85), and twelve types of eyepieces for visual observations and for taking photographs and motion pictures.

Using the electric circuitry of the picture-taking control console, the operator can do the following:

- turn on the lamp of the microscope illuminator;
- take motion pictures at a speed of 8 to 32 frames per second when using a type KCP-1M camera or at a speed of 8 to 64 frames per sec when using a type 16CP cine camera;
- take continuous sequences at speeds of 6 to 10 frames per sec;
- take sequences of single shots at intervals of 2 to 30 sec between the exposures;
- take sequences of single shots at intervals of 1.5 to 15 min between the exposure;
- take single shots.

BASIC CHARACTERISTICS


Magnification:

- in visual observations with binocular attachment 26 to $900\times$
- in visual observations with monocular attachment 17 to $600\times$
- in taking photographs 10.5 to $750\times$

in taking motion pictures with KCP-1M cine camera 6 to $412\times$

in taking motion pictures with 16CP cine camera 2.8 to $195\times$

Displacement range of mechanical stage in two
mutually perpendicular directions, mm 20×20
Accuracy of mechanical stage displacement readings, mm 0.1
Overall dimensions, mm $1400\times 750\times 1650$
Mass, kg 50



UNIVERSAL RESEARCH MICROSCOPE МБИ-15

The МБИ-15 biological microscope is a universal research microscope adapted for visual observation and taking photographs of objects by all existing microscopic methods.

With transmitted light, objects can be viewed against a bright field with direct and oblique illumination, against a dark field by the phase and interference contrast methods, and also with polarized light.


With reflected light, observations can be carried out on a bright and dark field, with mixed illumination, or using the fluorescence of the objects excited by the blue-violet region of the light spectrum in the range 400–450 nm and by ultra-violet rays up to 360 nm. Besides, the objects under observation can be illuminated from above and simultaneously from below using the phase contrast or dark-field technique. Photographs of objects are taken with a 35-mm film camera on 24×36 mm negatives and with a plate camera on 12×12 cm negatives. In taking photographs, the objects are illuminated from an impulse light source.

The microscope comes complete with seventeen objectives for transmitted light work, namely four planachromatic objectives (3.5×0.10; 9×0.20; 20×0.40; and 40×0.65), six apochromatic objectives (20×0.65; 40×0.95; 60×1.30; two 70×1.23; and 90×1.3), seven achromatic objectives (10×0.30; 10×0.40; 20×0.40; 30×0.90; 40×0.65; and two 90×1.25), three epi-objectives for reflected-light work (9×0.20; 21×0.40; and 40×0.65), and one achromatic objective (95×1.25). For visual observations and taking photographs, the microscope is furnished with a 4× Huygens eyepiece, a 7× Koelner eyepiece, compensating eyepieces (5×, 7×, 10×, 15×, and 20×) and a 12.5× orthoscopic objective.

The microscope is widely used in investigations in the fields of biology, botany, zoology and medicine, and also in other fields of science and technology.

BASIC CHARACTERISTICS

Magnification of microscope:		
in visual observations	• • • • •	28 to 2000×
in taking photographs	• • • • •	18 to 3000×
Overall dimensions, mm	• • • • •	1400×700×1140
Mass, kg	• • • • •	100



TRAVELLING BIOLOGICAL MICROSCOPE

"BIOLAM Д-3"

The "Biolam Д-3" travelling biological microscope is designed for the transmitted-light, bright-field examination of transparent objects in the fields of medicine, zoology, biology, and allied sciences under expeditionary conditions.

The optical system of the microscope has provisions for observation of objects through a type AY-12 binocular attachment.

The illuminating system is a demountable mirror which can be used for work with both natural (daylight) illumination and type OM-19 illuminator supplied on separate order.

The microscope comes complete with three achromatic objectives (20×0.40, 40×0.65, and 90×1.25), one planachromatic objective (9×0.20), and two paired compensating eyepieces (7× and 10×).

BASIC CHARACTERISTICS

Magnification of microscope 94.5 to 1350×
Focusing range, mm:		
coarse focusing 40
fine focusing 2
Overall dimensions, mm 230×140×350
Mass, kg 3.45



Stereoscopic Microscopes



STEREOSCOPIC MICROSCOPE MBC-2

The MBC-2 stereoscopic microscope on a universal stand forms an erect image of objects in solid relief which can be observed by both transmitted and reflected light. The microscope is used in the laboratories of medical, biological, zoological and other institutions in investigations on living objects, plants, etc.

The universal stand of the microscope is constructed so that it allows sufficient room for manipulation of objects which cannot be accommodated on the microscope stage.


The MBC-2 microscope is widely used in many other fields of science, technology and industry. Its objective is a special system built into an optical head and composed of four lenses with a focal distance of 80 mm and two paired Galileo systems. By rotating these systems, it is possible to obtain two ranges of magnification.

A distinction of the optical train of the MBC-2 microscope is a large working distance which remains the same for all magnifications, and a wide field of view. The objects observed in the MBC-2 microscope can be photographed with a type MΦH-5 photomicrographic attachment which is not included in the standard equipment of the microscope.

Observations with the MBC-2 microscope can be made with both natural (daylight) illumination and special-purpose illuminators.

BASIC CHARACTERISTICS

Magnification of microscope	3.5 to 88×
Linear field of view, mm	39 to 2.6
Working distance of microscope, mm	64
Magnification of optical head	0.57×, 1×, 2×, 4×, 7×
Magnification of eyepieces	6×, 8×, 12.5×
Magnification of eyepiece micrometer	8×
Overall dimensions, mm	400×245×460
Mass, kg	17.5



STEREOSCOPIC POLARIZING MICROSCOPE МПС-2

The МПС-2 stereoscopic polarizing microscope is designed for the investigation of petrographic sections of normal size (28×48 mm) and enlarged size (60×90 mm) by transmitted and reflected light (unpolarized and polarized). The microscope can also be used to investigate both polished and unpolished microsections with direct and oblique illumination. With the image produced by the МПС-2 microscope in solid relief, it is an easy matter to identify both transparent and some opaque minerals from their properties manifested in natural state, which fact considerably simplifies studies into the structural relationships between minerals. Alternating these observations with examinations of a polished microsection by reflected diffuse light, one can readily bring out most of the transparent and opaque minerals in a single section.

Owing to the long working distance of the microscope, a variety of manipulations can be done on the object.

Photographs of the image can be taken with a type МФН-5 photomicrographic attachment.

The microscope is used in mineralogy, petrography, mineralogy, crystallography, chemistry, physics, metallurgy, textile industry for fibre examinations, and some other fields of science and industry.

BASIC CHARACTERISTICS

Magnification of microscope	4 to 100×
Overall dimensions, mm	540×200×500
Mass, kg	20



BINOCULAR MICROSCOPE BM-51-2

The BM-51-2 binocular microscope is designed for examination of minute objects, engravings, surface patterns, fractures, breaks, scratches, and surfaces under small magnifications in cases where it is important to have a stereoscopic image.

Owing to the excellent stereoscopic effect produced by the microscope, details of objects can be viewed at any depth and of any shape.

The microscope is used in the assembly of precision mechanisms, and also in research work in medicine, botany, and zoology.

The BM-51-2 binocular microscope is in effect a combination of two self-contained optical trains each using a $0.7\times$ objective. The prisms of each optical train are built into a separate drawtube, and the top covers of the draw tubes mount eyepieces with a magnification of $12.5\times$.

BASIC CHARACTERISTICS

Total magnification of microscope	8.75×
Linear field of view, mm	25
Maximum working distance, mm	140
Overall dimensions, mm	dia. 160×300
Mass, kg	3.5



STEREOSCOPIC MICROSCOPE MBC-200


The MBC-200 is a stereoscopic microscope forms a direct image of the object in solid relief, by reflected light.

The illuminating system of the microscope produces shadowless illumination.

The microscope is widely used by the laboratories and shops of the microelectronics industry.

BASIC CHARACTERISTICS

Magnification of microscope	16.2 to 204×
Field of view, mm	7.3 to 1.0




STEREOSCOPIC MICROSCOPE MCCO

The MCCO stereoscopic microscope is intended for the examination of small objects, touching up photographic plates and other uses involving transmitted and reflected light, and mixed illumination.


The MCCO microscope can successfully be employed in the laboratories and shops of the microelectronics industry.

BASIC CHARACTERISTICS

Magnification of microscope	3.5 to 120×
Field of view, mm	44 to 1.8



Metallographic Microscopes



SIMPLIFIED METALLOGRAPHIC MICROSCOPE MMY-3


The MMY-3 simplified metallographic microscope is designed for visual examination of opaque objects by reflected light on a bright and a dark field, and by polarized light.

With the MMY-3 microscope, objects can be photographed with Type MΦH photomicrographic attachment.

The microscope is used at factory laboratories and for a variety of shop work.

BASIC CHARACTERISTICS

Magnification	100×, 300×, 500×
Range of stage displacements in two mutually perpendicular directions, mm	40×20
Maximum height of object, mm	100
Overall dimensions, mm	250×220×330
Mass, kg	6



40

METALLOGRAPHIC SHOP MICROSCOPE MMP-2


The MMP-2 is a metallographic shop microscope of inverted construction (that is, with the interchangeable stage located at the top), adapted to examine, inspect, and take photographs of the microstructure of metal and alloy microsections by the bright-field technique with direct and oblique illumination, and also with dark-field and polarized-light illumination. The optical systems of the microscope gives standard total magnifications necessary for visual observations with a binocular attachment, on a viewing screen, and for taking photographs of objects on 9×12 cm plates or on film with 24×36 mm negatives.

The automatically and manually actuated stages have motions in two mutually perpendicular directions. The objectives can be changed by rotating a revolving nosepiece.

The MMP-2 microscope is widely used in the shop laboratories of metal-working and metallurgical plants, at research institutions and other establishments concerned with the working and testing of metals.

BASIC CHARACTERISTICS

Magnification of microscope:	
in visual observation	80 to 1000×
in taking photographs	50 to 1000×
Displacement range of stage in two mutually perpendicular directions, mm	30×20
Accuracy of stage displacement readings, mm	0.1
Overall dimensions, mm	550×225×395
Mass, kg	24.5



41

METALLOGRAPHIC MICROSCOPE MMP-2P

The MMP-2P metallographic microscope is intended to examine and take photographs of the microstructures of metals, alloys and other opaque objects on a bright field by direct and oblique illumination, on a dark field, and with polarized light.

The optical system of the microscope provides standard magnifications necessary for visual observations with a binocular attachment, on a viewing screen, and in taking photographs of objects on 9×12 cm plates and on a film with 24×36 mm negatives.


In contrast to the MMP-2, the MMP-2P microscope comes complete with a stage manually controlled in two mutually perpendicular directions.

The objectives are changed with the aid of a revolving nosepiece.

The MMP-2P microscope is used in the shop laboratories of metal-working and metallurgical plants, at research institutions, and other establishments.

BASIC CHARACTERISTICS

Magnification of microscope:	
in visual observation	80 to 1000×
in taking photographs	50 to 1000×
Range of stage displacement in two mutually perpendicular directions, mm	
	30×20
Accuracy of stage displacement readings, mm	0.1
Overall dimensions, mm	550×225×395
Mass, kg	24.5



UNIVERSAL METALLOGRAPHIC MICROSCOPE МММ-9

The МММ-9 universal metallographic microscope is a research instrument designed to observe visually and to take photographs of the microstructures of opaque objects on a bright field with direct and oblique illumination, on a dark field with polarized light, by the phase contrast and interference method, to measure relief, and to evaluate fractures and microsections on a point-scale basis.

The optical system of the МММ-9 microscope incorporating planachromatic and planapochromatic objectives produces high-quality images of the objects.

Photographs are taken on 9×12 cm and 13×18 cm plates and on film with 24×36 mm negatives. The microscope has provisions for examination of objects on a viewing screen. The stage of the microscope has mechanisms for moving the object in a rectangular coordinate system.

The stage is lifted, the bellows of the camera are extended, and the exposure times are set from a control console.

The microscope is widely used in the metallographic laboratories of research institutions, in industry, and at educational establishments.

BASIC CHARACTERISTICS

Magnification of microscope	20×, 30×, 100–4000×
Range of stage displacement in two mutually perpendicular directions, mm	15×15
Size of viewing screen, mm	130×180
Overall dimensions, mm	1400×750×1200
Mass, kg	200



Polarizing Microscopes



RESEARCH POLARIZING MICROSCOPE МНН-8

The МНН-8 is a large polarizing microscope designed to examine transparent objects by transmitted unpolarized and polarized light in conoscopic and orthoscopic light beams.

The microscope comes complete with two planachromatic objectives (3.5×0.10 and 9×0.20), four achromatic objectives (20×0.40 ; 35×0.65 ; 60×0.85 , and 90×1.25), two sizes of Huygens eyepieces ($10\times$ and $8\times$), and two photographic eyepieces ($10\times$ and $15\times$).

The microscope is adapted to take a type ОИ-12 illuminator which may be used in the examination of opaque objects by reflected light (both polarized and unpolarized). When fitted with type МФН photographic attachments, the microscope can take photographs of the objects. Besides, the microscope can be fitted with a Fyodorov stage, a type КФ-4 attachment for phase-contrast work, and type И-13 dark-field condensor which are not however included in the standard equipment of the microscope and may be had on separate order.

The microscope is used by crystallographic, mineralogical, petrographic, chemical, physical, and many other research and educational laboratories.

BASIC CHARACTERISTICS

Magnification of microscope	17.5 to 1350 \times
Point displacement due to rotation of stage, mm	0-0.005
Range of stage rotation	360°
Accuracy of stage rotation readings	6'
Overall dimensions, mm	210 \times 390 \times 380
Mass, kg	8.4



TRAVELLING POLARIZING MICROSCOPE

МПД-1

The МПД-1 is a small polarizing microscope designed to examine mineralogical sections or powdered minerals and rocks by transmitted light under expeditionary conditions. Investigations can be carried out with both unpolarized and polarized light in orthoscopic and conoscopic light beams.

The microscope comes complete with a planachromatic objective (9×0.20), two achromatic objectives (40×0.65 and 60×0.85), and two sizes of eyepieces ($5\times$ and $8\times$).


The microscope is adapted for the investigation of minerals by the so-called immersion technique which is best suited to field conditions and reveals the optical characteristics of crystals sufficient to identify a given material.

Depending on the conditions at hand, examinations can be carried out with both natural (daylight) illumination and with special-purpose illuminators. The latter may be types ОИ-19 and ОИ-31 which are not included in the standard equipment of the microscope and may be had on separate order.

The microscope may take a type КФ-4 standard phase-contrast attachment and also a type ОИ-13 dark-field condensor.

BASIC CHARACTERISTICS

Magnification of microscope	45 to 900×
Overall dimensions, mm	290×200×100
Mass, kg	3.4



ORE MICROSCOPE МИН-9

The МИН-9 ore microscope is a polarizing microscope designed to examine opaque objects by polarized and unpolarized reflected light on a bright field.

It can also be used for examination of transparent objects by transmitted light at low magnifications with both unpolarized and polarized light.

The microscope comes complete with seven achromatic objectives (4.7×0.11; 9×0.20; 11×0.20; 21×0.40; 30×0.65; 40×0.65; and 95×1.25) and four different types of eyepieces (7×, 10×, 15×).

The microscope is mainly used by mineralographic and coal-petrographic laboratories, but can also be used by laboratories in other fields of science, technology and industry.

BASIC CHARACTERISTICS

Magnification of microscope	33 to 1425×
Least graduation of stage	1°
Accuracy of stage rotation readings	6'
Overall dimensions, mm	190×320×380
Mass, kg	5



Fluorescence Microscopes




FLUORESCENCE MICROSCOPE МЛ-2А

The МЛ-2А fluorescence microscope is intended to examine micro-biological, histological and other preparations by the fluorescence excited by the blue-violet region of the spectrum and also by rays with a wavelength of up to 360 nm.

In fluorescence examination, objects may be illuminated either from above through an opaque illuminator and the objectives or from below using a type ОИ-13 dark-field condensor or a type КФ-4 phase-contrast attachment included in the standard equipment of the instrument. The microscope may also be used to examine objects by mixed illumination. In this case, the light exciting fluorescence is directed simultaneously from above and from below; the preparation is illuminated by the phase-contrast or dark-field method.

The microscope may additionally be used for the examination of objects in the visible region of the spectrum by transmitted light on a bright field and by reflected light on a dark field.

The transmitted and reflected light objectives supplied with the microscope (10×0.40Л; 20×0.40Л; 30×0.90Л; 40×0.65Л; 40×0.75Л; 60×1.25Л; 60×1.0Л; 65×1.1Л; 90×1.25Л 21×0.40; 40×0.65; 40×0.30ФЛ; 40×0.65ФЛ; and 90×1.25ФЛ) produce a high quality image of the objects.



54

The МЛ-2А microscope comes complete with a type АУ-26 binocular attachment for examination of objects and a type МФН-10 monocular photographic attachment for taking photographs of objects with a ZORKY-4 photographic camera.

BASIC CHARACTERISTICS

Magnification of microscope:

in visual observations	53 to 1350×
in taking photographs	20 to 540×

Range of stage displacement in two mutually perpendicular directions, mm

40×60

Maximum height of object in reflected light work, mm


90

Overall dimensions, mm:

microscope	625×225×416
electric control panel	380×210×220

Mass, kg:

microscope	24
electric control panel	20.5



55

FLUORESCENCE MICROSCOPE МЛ-2Б


The МЛ-2Б fluorescence microscope is intended to examine micro-biological, hystological and other preparations by the fluorescence excited by the blue-violet region of the spectrum and also by ultra-violet rays with a wavelength of up to 360 nm.

In performance and design, the МЛ-2Б is almost identical with the МЛ-2А microscope.

The МЛ-2Б microscope comes complete with a trinocular attachment for observation and photomicrography and a type МФН-10 monocular photographic attachment to take photographs of objects with a ZORKY-4 photographic camera.

BASIC CHARACTERISTICS

Magnification of microscope:		
in visual observations	53 to 1350×	
in taking photographs	20 to 648×	
Range of stage displacement in two mutually perpendicular directions, mm		40×60
Maximum height of object in reflected light work, mm		90
Overall dimensions, mm:		
microscope	625×225×416	
electric control panel	380×210×220	
Mass, kg:		
microscope	24	
electric control panel	20.5	



FLUORESCENCE MICROSCOPE МЛ-2Б

The МЛ-2Б fluorescence microscope is intended to examine micro-biological, hystological and other preparations by the fluorescence excited by the blue-violet region of the spectrum and also by ultra-violet rays with a wave-length of up to 360 nm.


In performance and design, the МЛ-2Б is almost identical with the МЛ-2А microscope.

The МЛ-2Б has an additional facility for evaluating the intensity of fluorescence of the object. For this purpose, the standard equipment of the microscope includes a ФМЭЛ-1А photometric attachment with accessories.

The МЛ-2Б microscope is fitted with a type АУ-26 binocular attachment for observation of objects and a type МФН-10 monocular attachment to take photographs of objects with a ZORKY-4 photographic camera.

BASIC CHARACTERISTICS

Magnification of microscope:		
in visual observations	53 to 1350×	
in photomicrography	20 to 540×	
Range of stage displacements in two mutually perpendicular directions, mm		40×60
Maximum height of object in reflected light work, mm		90
Minimum area on preparation in photometry, mm		0.001
Overall dimensions, mm:		
microscope	625×225×416	
electric control panel	380×210×220	
Mass, kg:		
microscope	24	
electric control panel	20.5	



FLUORESCENCE MICROSCOPE МЛ-3

The МЛ-3 fluorescence microscope is designed to examine and take photographs of biological sections in the reflected light of fluorescence excited by the blue-violet region of the spectrum in the wave-length range 400–650 nm and by ultra-violet rays with a wave-length of up to 360 nm.


With the МЛ-3 microscope, objects can be examined by the light of fluorescence with illumination from above, through the opaque illuminator, and the objective by the bright field technique. In each case, the object can simultaneously be illuminated from below and examined by the dark field or phase contrast method, using a type ОИ-13 dark-field condensor or a type КФ-4 phase-contrast attachment. With the МЛ-3 microscope, photomicrographs can be taken on 35-mm films with 24×36 mm negatives using a photographic camera, photographic drawtube and an adapter supplied with the instrument. To measure the intensity of fluorescence, a type ФМЭЛ-1 photometric attachment may be fitted on the microscope head.

The type ОИ-13 condensor, КФ-4 phase contrast attachment, and the ФМЭЛ-1 photometric attachments are not included in the standard equipment of the microscope and may be had on separate order.

The МЛ-3 microscope is used by medical and biological laboratories, including hystological, cytological, microbiological, botanical, and many others.

BASIC CHARACTERISTICS

Magnification of microscope:		
in visual observation	• • • • •	50 to 1350×
in photomicrography	• • • • •	17 to 450×
Overall dimensions, mm	• • • • •	570×230×400
Mass, kg	• • • • •	30



TRAVELLING FLUORESCENCE MICROSCOPE МЛД-1

The МЛД-1 travelling fluorescence microscope is designed to examine visually biological preparations by the fluorescence excited by the blue-violet region of the spectrum in the wave-length range 400–550 nm and by ultra-violet rays with a wave-length of up to 365 nm.

With the МЛД-1 microscope, fluorescence examinations may be carried out only with objects illuminated from above, through a micro-objective, by the bright-field method.


For its operation, the microscope depends on the fluorescence given up by objects when excited by rays of a definite spectral composition.

The microscope comes complete with four achromatic objectives (10×0.40; 40×0.65 and 90×1.25) and four sizes of compensating eyepieces (5×, 7×, 10×, and 15×).

The microscope is mainly intended for field work in truck-mounted laboratories, but may also successfully be used in stationary laboratories concerned with medical and biological investigations.

BASIC CHARACTERISTICS

Excitation range fluorescence, Å	3600–4400
Fluorescence range investigated, Å	4000–6500
Magnification of microscope in visual observations	50 to 1350×
Range of stage displacement in horizontal plane, mm	60×40
Vertical displacement of stage, mm	40
Maximum height of object, mm	90
Overall dimensions of microscope in case, mm	570×230×400
Mass of microscope in case (less control panel), kg	27



CONTACT FLUORESCENCE MICROSCOPE МЛК-1


The МЛК-1 contact fluorescence microscope is designed to examine and take photographs of the microstructures of human and animal tissues by visible fluorescence in the wave-length range from 400 to 700 nm, excited by the blue-violet region of the spectrum. Among other things, the microscope can be used for the contact investigation of hard-to-reach human organs, for observation and photomicrography during operations and after removal of organs from the body and for diagnosing tumours and some other diseases. The object is illuminated with an opaque illuminator, from above, through the objective.

The microscope comes complete with a specially designed contact achromatic objective (11.6×0.65) and two specially designed eyepieces (8× and 15×). Photomicrographs are taken on 35-mm film with 24×36 mm negatives, using a strobe-lamp as a light source.

The МЛК-1 microscope is widely employed by medical research and clinical establishments in the field of surgery, oncology, gynecology, etc.

BASIC CHARACTERISTICS

Magnification of microscope	93× and 174×
Overall dimensions, mm:	
microscope	470×180×300
stand	800×650×1300
Mass, kg:	
microscope	5
stand	50



FLUORESCENCE MICROSCOPE МЛ-4

The МЛ-4 fluorescence microscope is designed to examine biological and other objects by visible and ultra-violet fluorescence. Fluorescence examinations may be carried out with the objects illuminated from above, through an opaque illuminator, and the objective.

Besides, the microscope may be used to examine objects illuminated from below, through a condensor, and also by the phase contrast technique.

The microscope has provisions for taking photographs of objects, for measuring the intensity of visible and ultra-violet fluorescence, and for registering the relative intensities of fluorescence in the various spectral intervals.

BASIC CHARACTERISTICS

Magnification of microscope:


in visual observations	50 to 1350×
in photomicrography	17 to 450×

Overall dimensions, mm:


microscope	730×500×550
control panel	320×400×900

Mass, kg:

microscope	50
control panel	72



Infra-red and Ultra-violet Microscopes



INFRA-RED MICROSCOPE MİK-4

The MİK-4 infra-red microscope is designed to investigate objects by unpolarized and polarized light in the infra-red region of the spectrum extending from 750 to 1200 nm in wave length, and also in the visible region of the spectrum.

Objects can be examined visually and their photographs can be taken both by transmitted and reflected light on a bright and dark field.

In transmitted-light work using polarized rays in the visible and infra-red region of the spectrum, use may be made of both orthoscopic and conoscopic light beams.


The infra-red image is converted to a visible one by a system incorporating an electron-optical image converter which has a magnification of 6 to 33.6 \times of its own.

The microscope comes complete with six objectives for transmitted-light work (3.5 \times 0.10; 9 \times 0.20; 60 \times 0.85; 10 \times 0.30; 75 \times 1.00; and 100 \times 0.80) and four objectives for reflected-light work (9 \times 0.20; 21 \times 0.40; 40 \times 0.65; and 95 \times 1.00). The standard equipment of the microscope also includes three types of projection eyepieces and five types of eyepieces (5 \times , 7 \times , 8 \times , 15 \times , and 28 \times) for visual observations.

The microscope may be used in mineralogy, study of semi-conductors, and also in glass-making.

BASIC CHARACTERISTICS

Magnification:	
in the infra-red region	44 to 5517 \times
in the visible range	21 to 1425 \times
Range of object displacement, mm	38 \times 22
Overall dimensions, mm	500 \times 175 \times 375
Mass, kg	17.5



ULTRA-VIOLET MICROSCOPE MYΦ-6

The MYΦ-6 ultra-violet microscope is designed to investigate and make photographs of biological preparations by transmitted light and to produce photographic records of optical density and absorption spectra of the various areas of preparations in the region from 2500 to 7000 Å.

The microscope also makes it possible to study biological preparations by their visible and ultra-violet fluorescence. Ultra-violet fluorescence is investigated by photomicrography and subsequent application of the colour transformation method with the aid of a type MYΦ-2X chromoscope which is part of the standard equipment of the MYΦ-6 microscope.

The microscope comes complete with two ultra-violet objectives (10×0.20 and 58×0.80), one quartz objective (48.5×0.65); five mirrors objectives (40×0.50; 75×0.65; 73.5×1.0; 115×0.7; 125×1.0) and four achromatic objectives (10×0.40; 30×0.90; 40×0.75, and 60×1.25).

The microscope is also supplied with four Huygens eyepieces (5×, 7×, 10×, and 15×) and nine different projection quartz eyepieces (3×, 8×, 10×) for photomicrography. Visual observation of preparations by transmitted ultra-violet rays is carried out on the screen of an electron-optical image converter.

The MYΦ-6 microscope is widely used by various medical and educational establishments, and also by research institutions in the fields of biology, botany, virusology and similar divisions of science.

BASIC CHARACTERISTICS


Magnification in visual observation:

by ultra-violet rays on the screen	
of an electron-optical image converter	28 to 948×
by visible light	40-1350×
by visible fluorescence	40 to 630×

Magnification in photomicrography:

by ultra-violet rays	28 to 948×
by visible fluorescence	28 to 853×
by ultra-violet luminescence	28 to 440×

Overall dimensions, mm	1230×350×650
Mass, kg	200



INFRA-RED STEREOSCOPIC MICROSCOPE

ММК-3

The ММК-3 infra-red stereoscopic microscope is designed for stereoscopic examination of objects transparent in the wave-length range from 0.75 to 1.2 μ of the infra-red region of the spectrum.

The microscope has provisions for visual stereoscopic and non-stereoscopic observations of objects and for taking their photographs in infra-red light.

Photographs are taken with a ZORKY-4 camera on a 35-mm film with 24 \times 36 mm negatives.

The infra-red image is converted to a visible one by electron-optical image converters.

The microscope is employed in the investigations of semi-conductors, in glass-making, mineralogy, petrography, zoology, etc.

BASIC CHARACTERISTICS


Magnification of microscope:

in visual observations 25 to 544 \times


in photomicrography 6 to 120 \times

Overall dimensions, mm 470 \times 180 \times 450

Mass, kg 15



Microprojection, Photomicrography and Cine-micrography Outfits



TRICHINA MICROPROJECTOR ТМП

The ТМП trichina microprojector is designed to examine and detect trichina and other parasites in meat products.

The trichonoscope comes complete with a 30× achromatic objective.

The preparation to be examined is squeezed between two special plates directly on the microscope.

The stage movement makes it possible to investigate all rows of the squeezer cells in turn. A special device prevents examination of cells in the second row before all cells in the first row have been examined.

The instrument has a system of limit stops, owing to which the squeezer may be removed only in two positions, namely either before the beginning of a test or after all the cells of the squeezer have been examined.

The microprojector is widely used by the laboratory of meat-packing plants, meat inspectors on markets, etc.

BASIC CHARACTERISTICS

Size of image on screen (ellipse with axes), mm	280 to 300
Distance from screen to observer, mm	350
Displacement range of stage and squeezer, mm:	
in longitudinal direction	150
in transverse direction	32
Vertical travel range of objective (change of focus), mm	±2
Overall dimensions, of packing case, mm	1780×1050×815
Mass of microprojector, kg	70




PHOTO MACRO- AND MICROGRAPHY OUTFIT ФМН-2

The ФМН-2 outfit is designed to take photographs of various transparent and opaque objects encountered in research work in the fields of technology, biology, botany, etc.

The outfit has also provisions for making reproductions of printed matter, hand drawings, mechanical drawings, and the like.


In conjunction with biological microscopes, the ФМН-2 outfit may be used to take photomicrographs of transparent objects. The drawing attachment supplied with the outfit (model ФМН-2Р) makes it possible to draw microobjects on paper.

The standard equipment of the ФМН-2 outfit includes special-purpose illuminators which can produce direct, oblique, shadowless, annular, and all other forms of illumination for both transparent and opaque objects. With these illuminators, any object may be illuminated with diffuse or condensed light.

The ФМН-2 outfit comes complete with a set of photographic objectives with focal distances of 40, 65, 100 and 150 mm.

BASIC CHARACTERISTICS

Size of photographic plates, cm	9×12
Range of scales in photography	0.5 to 20×
Maximum size of object, mm:	
in operation with small macrostage	50 dia
in operation with large macrostage	100 dia
Overall dimensions, mm	1200×550×600
Mass, kg	60



PHOTOMICROGRAPHIC OUTFIT ФМН-3

The ФМН-3 outfit is designed to take photographs of objects through a microscope with the use of a high-intensity light source.

Photographs are taken on photographic plates measuring 9×12 cm or on 35-mm film with 24×36 mm negatives.

The light filters supplied with the outfit are placed in pockets on the outside of the field-diaphragm assembly. Any other light filters available to the investigator can be installed in a pocket on the heat-filter mount. The microscope is set up on the base of the ФМН-3 outfit and is locked in place by angles and screws.

The outfit is intended for use at research and educational laboratories in conjunction with biological microscopes.

BASIC CHARACTERISTICS

Size of photographic plate, cm	6.5×9 and 9×12
Size of negative on film, mm	24×36
Overall dimensions, mm	300×300×800
Mass, kg	20



76

CINEMICROGRAPHY OUTFIT МКУ-2


The МКУ-2 is a stationary universal outfit incorporating facilities for the examination of microscopic objects by the various light-microscopy methods and cinemicrography of fast and slow processes in microscopic objects over a broad range of shooting speeds.

The outfit will be of interest to biological, medical, crystallographic and many other laboratories at research establishments, educational institutions and industrial enterprises.

BASIC CHARACTERISTICS



Magnification of the optical system (without microscope) in cinemicrography	0.76 to $16 \times$
Continuous shooting speed of fast processes	10 frames/s to 64 frames/s

Single-shot micrography of slow processes at intervals of	from 2 s. to 15 min.
Overall dimensions, mm	700×1300×1800
Mass, kg	110



77

**Microscopes
for Advanced Microscopic
Investigations**



MICROSPECTRAL ANALYSIS OUTFIT MCJI-2

The MCJI-2 outfit is designed for the microspectral emission local analysis of structural components of alloy metals, and minerals.

The outfit comprises a microscope, an electrode holder, a laser, spectrograph, and a control console.


The microscope is used to examine objects to locate the areas to be analysed, and to project the laser emission onto the area under investigation. The microscope can be used to examine transparent and opaque objects on a dark field and with polarized light.

The laser serves to evaporate the material within the area under investigation.

The evaporated material is excited by a spark discharge.

The spectrograph decomposes the emission of the material under investigation into a spectrum and takes photographs of the latter.


The control console consists of a power unit and a system to select the various modes of operation for the laser and spark generator.



80

BASIC CHARACTERISTICS

Magnification of microscope in visual observation	50 to 400×
Range of stage displacement in two mutually perpendicular directions, mm	40×20
Wave-length of laser emission, μ	1.06
Effective range of spectrograph, μ	0.22 to 0.88
Power taken by control console, W	800
Least diameter of area under analysis, μ	40
Overall dimensions, mm:	
microscope	560×240×540
laser	390×130×130
control console	545×570×1550
spectrograph	1000×710×400
desk	1860×700×760
Mass, kg	424



81

CRIMINAL INVESTIGATION COMPARISON MICROSCOPE MCK-1

The MCK-1 criminal investigation comparison microscope provides a means for the examination and photography of two objects under comparison which can be seen simultaneously or separately in the field of vision.

The optical system of the microscope consists of two identical and independent optical trains forming an image of two objects in the same plane. In observations with a binocular attachment, projection on a viewing screen, and photography, the microscope forms erect images of two objects under comparison within the same field. The magnification can be selected by bringing in or out a particular Galileo system which gives four magnifications as follows: 2.5, 0.6, 0.4, and 1.6 \times .

Using the light filters supplied with the microscope, objects can be examined by fluorescence.

Photographs are taken on 6.5 \times 9 cm plates.

The construction of the MCK-1 microscope ensures an ultimate in convenience to the observer and reliability of operation.

The microscope is widely used in criminal investigation laboratories.

BASIC CHARACTERISTICS

Magnification of microscope:


in visual observation	2.7 to 60 \times
in photography	2 to 21 \times

Size of viewing screen, mm 130 to 150

Displacement range of stage, mm 40 \times 40

Overall dimensions, mm 800 \times 500 \times 850

Mass, kg 80



INTERFERENCE MICROSCOPE МБИН-4

The МБИН-4 interference biological microscope is designed to examine objects by transmitted light using the interference contrast method, to measure the path length difference due to films and sections, and to determine from this path length difference the dry weight of living and fixed cells and sections, dry matter concentration, the total weight, thickness and volume of an object.

The microscope comes complete with two objectives (10×0.30 and 40×0.65) and four sizes of eyepieces ($7\times$, $10\times$, $15\times$ and $20\times$).

In visual observation and photomicrography with the aid of type МФН-11 binocular attachment, the magnification may be adjusted by rotation of a revolving nosepiece by $1.1\times$, $1.6\times$, $2.5\times$.

The microscope is employed at factories laboratories, research institutions, and educational establishments.

BASIC CHARACTERISTICS

Magnification of microscope:


in visual observation	77 to 960 \times
in photomicrography	26 to 240 \times

Size of split image, μ	165 and 330
----------------------------	-------------


Sensitivity of interference method, λ	0.02 to 0.003
---	---------------

Overall dimensions, mm	600 \times 150 \times 500
------------------------	-------------------------------

Mass, kg	20
----------	----



Measuring Microscopes



SMALL TOOL-ROOM MICROSCOPE ММН-2

The ММН-2 small tool-room microscope is designed to measure linear dimensions in rectangular and polar coordinates, and also angle elements in various parts, cutting and measuring tools.


The ММН-2 tool-room microscope has provisions for measuring the elements of external threads, angles, tapers, radii, and effective dimensions of various templates, gages, and similar objects.

The microscope has provisions for observation with both transmitted and reflected light. Owing to a prism, an erect (non-inverted) image of the object under observation can be seen in the eyepiece, and all movement of the object are seen by the observer as they actually occur, which adds to the operator's skills.

The ММН-2 microscope is widely used in the shops and measurement laboratories of engineering, instrument-making and tool-making factories, at educational and research establishments. In shops, the ММН-2 microscope is the main measuring instrument in the manufacture of parts.

BASIC CHARACTERISTICS

Magnification of microscope	10×, 30×, 50×
Ranges of measurement, mm:	
in longitudinal direction	0-75
in transverse direction	0-25
Limit of error in length and diameter measurements, mm	±0.005
Maximum separation between centres, mm	200
Maximum diameter of objects mounted on centres, mm	55
Maximum diameter of objects set up on knife-edge supports, mm	100
Overall dimensions, mm	375×350×430
Mass, kg	27



LARGE TOOL-ROOM MICROSCOPE БМИ-1

The БМИ-1 large tool-room microscope is intended to measure the outside linear dimensions and diameters of shafts, the angles of cutting tools and templates, the dimensions of thread cutting tools and gauges, centre to centre distances of holes, dimensions of cylindrical and taper bushings, radial shapes, and similar components.

The microscope can also be used for marking out with a centre-punch.

With the БМИ-1 microscope, measurements may be made both directly with micrometers and by comparison when a magnified image of the outline to be measured is projected and compared with a reference outline on a drawing.


Threads can be measured by an outline comparison method, the axial section method, and with measuring knife edges.

The microscope has provisions for measurements in rectangular and polar coordinates, with transmitted and reflected light.

The multitude of measurements that can be made with БМИ-1 tool-room microscope makes it an indispensable instrument for tool-making shops, measuring laboratories, research establishments, and educational institutions.

BASIC CHARACTERISTICS

Magnification of microscope	10×, 15×, 30×, 50×
Ranges of measurement, mm:	
in longitudinal direction	0-150
in transverse direction	0-50
Maximum diameter of object mounted on centres, mm	85
Maximum diameter of part set up on knife-edge supports, mm	130
Maximum distance between objective and workable, mm	200
Overall dimensions, mm	840×800×870
Mass of microscope complete with angle-measuring eyepiece head and illuminator, kg	65




UNIVERSAL MEASURING MICROSCOPE YMM-23

The YMM-23 universal measuring microscope is a sighting instrument with projection-type reading devices designed to measure the linear and angle dimensions of various parts in rectangular and polar coordinates. Among other things, it can be used to measure all kinds of workpieces, cutting tools, templates, cams, tapers, threading taps and dies, chasers, etc. It can also be used to measure the dimensions of small holes.

The YMM-23 universal measuring microscope is widely used at engineering works, research institutions, and educational establishments.

BASIC CHARACTERISTICS

Range of length measurement, mm:	
in longitudinal direction	0-200
in transverse direction	0-100
Range of angle measurements	0-360°
Range of whole diameter measurement by contactless methods, mm	0.1-50
Overall dimensions, mm	1150×1150×1500
Mass, kg	541



UNIVERSAL MEASURING MICROSCOPE YMM-24


The YMM-24 universal measuring microscope is designed to measure the linear and angular dimensions of various parts in rectangular and polar coordinates. Among other things, can be used to measure all elements of thread gauges, templates, cutting tools, milling cutters, and similar parts.

The YMM-24 microscope is a measuring sighting instrument with projection-type reading devices. The workpiece being measured is held stationary, and the sighting system of the instrument together with the carriages is free to be moved in a horizontal plane in two mutually perpendicular directions.

The universal measuring microscope can be widely used in the laboratories of medium and heavy engineering works, and also at research establishments.

BASIC CHARACTERISTICS

Range of length measurement, mm:	
in longitudinal direction	0-500
in transverse direction	0-200
Range of angle measurements	0-360°
Maximum distance between centres, mm	1000
Maximum mass of workpiece to be measured, kg	100
Range of travel of lifting device, mm	250
Overall dimensions, mm	1500×1300×1600
Mass, kg	1550



READING MICROSCOPE МИР-2

The МИР-2 reading microscope is a simplified modification of a measuring microscope with a reading eyepiece scale and is used to measure small objects and distances between mark lines, dots, and surface irregularities.

The МИР-2 microscope is mainly used by shop and factory laboratories, quality inspection departments, and educational establishments.

BASIC CHARACTERISTICS

Magnification of microscope	19 to 33×
Range of measurement, mm	0.015 to 6
Overall dimensions, mm	134×34×67
Mass, kg	0.29



READING MICROSCOPE МПБ-2

The МПБ-2 reading microscope is used to measure the diameter of the indentation produced on metals in the Brinell hardness test.

The microscope is supplied with conversion charts from which the respective Brinell hardness number can be found for the diameter of the indentation measured with the microscope.

BASIC CHARACTERISTICS

Magnification of microscope	24×
Linear field of view, mm	9
Range of eyepiece adjustment, diopters	±4
Maximum diameter of indentation measured on scale, mm	6.5
Overall dimensions of packing case, mm	64×75×220
Mass of instrument, in packing case, kg	0.8



HARDNESS TEST MICROSCOPE МПВ-1

The МПВ-1 hardness test microscope is an accessory to hardness testers and is used to measure the linear dimensions of indentations in the Vickers hardness test.

The set of objectives and eyepiece supplied with the microscope provide for visual observations with magnifications of $50\times$ and $125\times$.

BASIC CHARACTERISTICS

Magnification of microscope	$50\times$ and $125\times$
Range of measurement, mm	0-12
Accuracy of reading, mm:	
with 10×0.30 objective	0.001
with 3.7×0.11 objectives	0.0025
Overall dimensions, mm	$210\times 165\times 95$
Mass, kg	2.29



94

MICROSCOPE FOR SPECTROGRAM MEASUREMENTS МИР-12

The МИР-12 microscope is used to interpret spectrograms by measuring distances between the unknown and known spectral lines.

In identifying the spectral lines of an unknown element, the wavelength "yardstick" to be used may be the spectrum of pure iron because the spectral lines of iron have been investigated in detail and their wave-lengths can be taken from spectral tables. For this purpose it is necessary to take photographs a spectrum of pure iron and juxtapose it edge to edge with that of the element under investigation.

The instrument is used in all laboratories and establishments concerned with spectral analysis.


BASIC CHARACTERISTICS

Magnification of microscope	$15\times$
Linear field of view, mm	5
Range of measurements in longitudinal direction, mm	0 to 50
Overall dimensions, mm	$440\times 280\times 315$
Mass, kg	12



95

**Surface
Inspection
Instruments**




MICROINTERFEROMETER ММН-10

The ММН-10 microinterferometer is used to measure the height of irregularities on machined surfaces.

With the ММН-10 microscope, measurements can be made by the conventional interference method ($0.65-0.03 \mu$) and also by the immersion-replica method ($10-0.6 \mu$).

The interference method is usually employed to measure reflecting machined surfaces with white and monochromatic light. The immersion-replica method is employed only in conjunction with monochromatic light to measure the height of irregularities on machined metallic and nonmetallic surfaces widely differing in reflectance, also to measure the height of irregularities on hard-to-reach or internal surfaces on parts of any size and shape without removal from a machine-tool.

The interference pattern is photographed with a ZORKY-4 camera on negative measuring 24×36 mm.



BASIC CHARACTERISTICS

Range of measurement of irregularities, μ 10-0.03

Magnifications:


in visual observations with 20× eyepiece micrometer	500×
in visual observations with 15× eyepiece	375×
in photography	200×

Field of view, mm:

in visual observations with 20× eyepiece micrometer	0.36
in visual observations with 15× eyepiece	0.40
in photography	0.12×0.18

Overall dimensions, mm 300×330×380

Mass, kg 30




MULTIBEAM MICROINTERFEROMETER ММН-11

The ММН-11 multibeam microinterferometer is a reflected-light microscope in the object plane of which is placed an interference device. The ММН-11 microinterferometer is used for the visual evaluation and measurement of the height of steps, thin film thickness, small linear displacements, and minute deviations from an ideal plane from the value of shift or curvature of interference bands in the range from 1 to 0.005 μ . The ММН-11 microinterferometer can measure surfaces with a reflectance of at least 0.3.


To ensure an optimum contrast of interference bands on surfaces widely differing in reflectance, the microinterferometer comes complete with four plates having reflectances of 0.45, 0.65, 0.75 and 0.90.

The microinterferometer may be used at the laboratories of research institutions and educational establishments, and also by industrial laboratories.



BASIC CHARACTERISTICS

Range of measurement, μ	1-0.005
Magnification of microinterferometer	50×
Field of view, mm	3.3
Overall dimensions, mm	170×350×380
Mass, kg	20




INTERFERENCE MICROSCOPE ММН-4

The ММН-4 interference microscope is a special-purpose high-precision measuring instrument used for the visual evaluation, measurement and photography of irregularities on the external surfaces of parts with parallel machining marks, and also the height of steps, film thickness, coating thickness, and similar objects.

The interference microscope is widely used at the laboratories of research institutions and educational establishments and industrial organization concerned with surface finish problems.

BASIC CHARACTERISTICS

Magnification of microscope:	
in visual observation	490×
in photography	290×
Field of view of microscope, mm:	
in visual observation	0.32
in photography	0.10
Range of measurements in surface finish control, μ	1.0 to 0.03
Overall dimensions, mm	340×300×380
Mass, kg	23.4



MICROPROFILOMETER ММН-12


The ММН-12 microprofilometer is a special-purpose instrument designed for visual evaluation and measurement of the height of irregularities in the range from 0.03 to 1 μ on outside surfaces on which machining marks occur at random.


The ММН-12 microprofilometer is an interference-type instrument in which surface measurements are made by using bands of the same chromatic order ("a fluted spectrum").

The instrument may be used at industrial laboratories, research institutions for the control of metallic and other well reflecting surfaces worked by honing or electrical polishing, and also to test surfaces on which machine marks are of local nature, such as after polishing of glass, quartz, germanium, silicon and similar materials with free abrasives.

BASIC CHARACTERISTICS

Range of measurement of surface regularities, μ	0.03-1
Magnification of the instrument	510×
Effective length in the plane of object, mm	0.25
Objective:	
focal distance, mm	10
numerical aperture	0.5
Overall dimensions, mm	250×500×380
Mass, kg	30





Внешторгиздат. Изд. № 253213